amieamed®








EXCEED



Gebrauchsanweisung (DE)
Operating manual (EN)
Gebruiksaanwijzing (NL)
Instrucciones de uso (ES)
Manual de instruções (PT)
Opskrift (DA)
Bruksanvisning (SV)
Instrukce (CS)
Инструкции за употреба (BG)

تعليمات الاستخدام (AR)
Kullanım kılavuzu (TR)
คู่มือการใช้งานฉบับตั้ง (TH)
Instrukcja obsługi (PL)
Notice d'utilisation (FR)
Notificare de utilizare (RO)
Notice d'utilisation (IT)
Käyttöilmoitus (FI)

CE 0123

	1	2*	3a	3b	4	5*	6	7	8
DE	Steuergerät	Hygienemodul	Handstück	Anschlusskabel	Handstückablage	Hygieneschlauch	Fußschalter	Stecker-netzteil	Länderadapter
EN	Controller	Safety cartridge	Handpiece	Handpiece cable	Handpiece tray	Handpiece cover	Foot switch	Power supply	Country-specific adapter
NL	Stuurapparaat	Hygiënemodule	Handstuk	Aansluitkabel	Handstukhouder	Hygiëneslang	Voetschakelaar	Stekkeradapter	Landadapter
ES	Dispositivo de mando	Módulo de higiene	Unidad de mano	Cable de conexión	Bandeja de unidad de mano	Funda higiénica	Interruptor de pedal	Adaptador de red	Adaptador de país
PT	Aparelho de controle	Módulo de higiene	Pen	Cabo de conexão	Base da pen	Manga de higiene	Interruptor de pedal	Transformador	Adaptador dos países
DA	Styreenhed	Hygiejнемодул	Håndsæt	Tilslutningskabel	Håndsætholder	Hygiejneslange	Fodkontakt	Strømforsyningsstik	Landespecifik adapter
SV	Styrenhet	Hygienmodul	Handenhet	Ström-kabel	Handenhetshållare	Hygieneslang	Fotbrytare	Kontakt-nättdel	Landsadapter
CZ	Řídicí jednotka	Hygienický modul	Ruční přístroj	Připojovací kabel	Uložení ručního přístroje	Hygienická hadice	Nožní pedál	Sifový adaptér	Adaptéru dané země
BG	Контролер	Хигиенен модул	Накрайник	Свързващ кабел	Поставка за крайник	Предпазен ръкав	Крачен превключвател	Мрежов захранващ блок	Специфичен за страната адаптер
AR	وحدة تحكم	وحدة النظافة	قلم الوخز	كابل التوصيل	حامل قلم الوخز	خرطوم نظافة	مفتاح تشغيل بالقدم	وحدة قدرة قابسة واحدة	محول كهربائي دولي واحد لكل منها
TR	Kumanda cihazı	Hijyen modülü	El parçası	Bağlantı kablosu	El parçası rafı	Hijyen hortumu	Ayak şalteri	Fiş adaptörü	Ülke adaptörü
TH	ชุดควบคุม	โมดูลสุขอนามัย	ด้ามจับ	สายเชื่อมต่อ	ที่วางด้ามจับ	ที่ครอบด้ามจับ	สวิตช์เท้า	อุปกรณ์จ่ายไฟ	อะแดปเตอร์เฉพาะสำหรับประเทศ
PL	jednostka sterująca / sterownik	moduł higieniczny	rekojeść	kabel przyłączeniowy rekojeści	uchwyt / podstawka na rekojeść	waż higieniczny	przetacznik nożny / pedał	zasilacz wtyczkowy	adapter krajowy (do gniazdek z różnymi standardami)
FR	unité de commande	module d'hygiène	pièce à main	câble de connexion de la pièce à main	support pour la pièce à main	tuyau d'hygiène	interrupteur à pied / pédale	bloc d'alimentation enfichable	adaptateur de voyage / adaptateur secteur international
RO	unitate de control / controler	modul de igienă	mâner	cablu de conectare a mânerului	mâner / suport pentru mâner	furtun de igienă	comutator de picior / pedală	alimentator plug-in	adaptor de țară (pentru prize de cu standarde diferite)
IT	unità di controllo / dispositivo di comando	modulo igienico	manipolo	cavo di collegamento del manipolo	supporto / base per il manipolo	tubo igienico flessibile	interruttore a pedale / pedale	alimentatore a spina	adattatore nazionale (per prese con standard diversi)
FI	ohjauksyksikö / ohjain	hygieniamoduuli	kahva	Kahvan liitäntäkappeli	Kahvan pidike / jalusta	hygienialetku	jalkakytin / poljin	virtalähde pistoke	maakohtainen sovitin (eri standardien mukaisille pistorasioille)
Artikelnummer	AD61MP	EMN-PAE-D06MP	CME21MP	5E-G770	GAD60	E-0610	E1038MP	E116503	EU: E-1154-W UK: E-1155-W AU: E-1157-W
Article description		(6-MNP 1.5-Exceed)							
Artikelnummer									
Número de artículo									
Referência de artigo									
Varenummer									
Artikelnummer									
Výrobní číslo									
Номер на артикул									
سم السلعة									
Ürün numarası									
หมายเลขสินค้า									
Numer katalogowy									
Référence									
Număr de catalog									
Numero di catalogo									
Katalogin numero									
		* Nicht im Lieferumfang enthalten * Not included with the delivery * Niet meegeleverd * No incluido en el alcance del envío * Não incluído no escopo de fornecimento * Medfølger ikke * Ingår inte i leveransinnehållet * Dodávka neobsahuje * Не се съдържа в обема на доставката							
						* غير مدرج في الملحقات الأساسية * Teslimat kapsamına dahil değildir * ไม่อยู่ในขอบเขตการส่งมอบ * Nie dołączone do zestawu * Non compris dans la livraison * Nu este inclus în set * Non incluso nel set * Ei sisälly sarjaan			

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CONTENTS

- 1 About this manual 32**
 - 1.1 Presentation of warning symbols 32

- 2 Important safety instructions 32**
 - 2.1 General safety instructions 32
 - 2.2 Product-specific safety instructions 32
 - 2.3 Important health and safety precautions 33
 - 2.4 Contraindications and side effects 34
 - 2.5 User requirements (qualifications)..... 35
 - 2.6 Intended purpose, indications and proper use 35
 - 2.7 Symbols on the product 35

- 3 Scope of delivery 37**

- 4 Product information..... 37**
 - 4.1 Technical data 37
 - 4.2 Operating conditions 37
 - 4.3 Accessories and spare parts..... 37

- 5 Getting started 38**
 - 5.1 Positioning the device 39
 - 5.2 Connecting the power supply 39
 - 5.3 Connecting the foot switch (optional)..... 40
 - 5.4 Connecting the handpiece 40
 - 5.5 Cleaning and disinfecting the equipment..... 40
 - 5.6 Pulling the handpiece cover on 41
 - 5.7 Inserting or replacing the safety cartridge..... 41
 - 5.8 Checking equipment..... 42

- 6 Use 42**
 - 6.1 Switching the control unit to operating mode or standby mode..... 42
 - 6.2 Configuring the controller (menu) 42
 - 6.3 Switching the handpiece on and off..... 43
 - 6.4 Setting the penetration frequency..... 44
 - 6.5 Setting the needle protrusion depth..... 44

- 7 Cleaning and maintenance..... 45**
 - 7.1 Routine maintenance..... 45
 - 7.2 Material compatibility..... 45
 - 7.3 Disinfecting surfaces..... 45
 - 7.4 Cleaning surfaces 46

8	Transport and storage conditions	46
8.1	Controller, handpiece with handpiece cable, power supply and foot switch	46
8.2	Safety cartridge	46
9	Disposal	46
10	Questions and problems	47
11	Manufacturer declarations	48
11.1	Guarantee.....	48
11.2	Electromagnetic compatibility.....	48

1 ABOUT THIS MANUAL

This operating manual applies to the EXCEED (AD61MP) and its accessories. It contains important information about setting up, operating and looking after the device safely and as intended for use.

The operating manual does not contain all information necessary for safe operation of the device and its accessories. Therefore, the following documents also need to be taken into account:

- Safety datasheets for disinfectants and cleaning materials
- Safety at the workplace requirements and statutory provisions for microneedling

1.1 Presentation of warning symbols

Warning symbols highlight the danger to persons or objects and are set up as follows:

 SIGNAL WORD	
Type of danger	
Consequences	
▶ Prevention	
Signal word	Indicates the seriousness of the hazard (see following table)

Type of hazard the type and source of the hazard

Consequences	describes possible consequences if ignored
--------------	--

Prevention	indicates how the hazard can be prevented
------------	---

SIGNAL WORD	MEANING
-------------	---------

Danger	Indicates a danger which will certainly result in severe or fatal injuries if the danger is not avoided
--------	---

Warning	Indicates a danger which could result in severe or fatal injuries if the danger is not avoided
---------	--

Caution	Indicates a danger which may cause slight to medium-severe injuries if the danger is not avoided
---------	--

Attention	Indicates a potential risk of damage to the environment, property or the equipment if the danger is not avoided
-----------	---

Symbols in this operating manual

SYMBOL	MEANING
▶	Action required
•	Bullet points
-	Sub-bullet points

2 IMPORTANT SAFETY INSTRUCTIONS

2.1 General safety instructions

- ▶ Read this operating manual carefully and in full.
- ▶ Keep this operating manual in a location which is accessible to all those using, cleaning, disinfecting, storing or transporting the device.
- ▶ Always include this operating manual when transferring ownership of the device.
- ▶ Follow the safety regulations for the microneedling in your country. Keep your workplace hygienically clean and ensure that you have adequate lighting.
- ▶ Only use the device, its accessories and all handpiece cables when in perfect condition.

- ▶ Only use original safety cartridges, accessories and spare parts from amiea med.

2.2 Product-specific safety instructions

- ▶ Never modify the device, the safety cartridges or other accessories.

Modifications can lead to an electric shock or a needle puncture wound with a contaminated needle.

- ▶ Prevent fluids from getting inside the device, the handpiece, the foot switch or the power supply.
- ▶ Position the device such that the recommended protective distance from possible sources of interference is maintained. You can find information on the recommen-

ded protective distances in the chapter Recommended protective distances from portable and mobile HF communications devices. Avoid sources of interference by ensuring that portable or mobile radio devices are not operated in the vicinity of the device.

- ▶ When handling the handpiece and the handpiece cable, protect them from contamination by bodily fluids or substances contaminated with bodily fluids by means of a handpiece cover (see chapter "5.7 Inserting or replacing the safety cartridge" on page 41).
- ▶ If the device is not being used, disconnect it from the power supply and place the handpiece in the handpiece tray so that it cannot roll away or fall.
- ▶ Observe the technical data provided in this operating manual and comply with the operating, transport and storage conditions (see chapter "8 Transport and storage conditions" on page 46).
- ▶ Return the device to a specialist retailer for inspection if there are signs of damage, if it does not function as normal or if liquids enter the device or handpiece.
- ▶ Return the device to a specialist retailer for inspection in accordance with the manufacturer's instructions (safety check), at regular 24 month intervals. In this regard, refer to chapter Inspection.
- ▶ Before use, always ensure that the date for the next inspection (safety check) on the test plate of the control unit has not passed. Do not use the device if the date on the test plate has passed and commission a new inspection if necessary.

2.3 Important health and safety precautions

Please observe the following instructions in order to prevent contaminations or infectious diseases being transmitted to the patients or users during treatment:

- ▶ Follow all of the steps for disinfecting the equipment before treatment (see chapter "5.5 Cleaning and disinfecting the equipment" on page 40).
- ▶ Use disposable gloves made from nitrile or latex during the treatment and disinfect these before use. Observe the applicable directives for your country when selecting suitable disinfectants.
- ▶ Clean the respective area of the patient's skin with a mild cleaning agent and disinfectant before treatment. Observe the applicable directives for your country when selecting suitable disinfectants.
- ▶ During use, the handpiece, the handpiece cable and the device including the integrated handpiece tray must be completely covered with protective film. The handpiece

must be covered with the protective film before attaching the safety cartridge (see chapter "5.6 Pulling the handpiece cover on" on page 41).

- ▶ safety cartridges are sterile disposable products (consumables) and are only ever permitted to be used once!
- ▶ Dispose of used or faulty safety cartridges as well as safety cartridges whose packaging is damaged in a non-penetrable sharps container in accordance with the regulations applicable in your country.
- ▶ Before the treatment, switch the handpiece off and check that the needles have been completely withdrawn into the safety cartridge. If this is not the case, dispose of the safety cartridge immediately.
- ▶ Ensure that safety cartridges never come into contact with contaminated items, such as clothing for example. Contaminated safety cartridges must not be used but rather disposed of immediately.
- ▶ Check at regular intervals whether the handpiece is visibly contaminated by bodily fluids running back or by the accompanying preparations being used. In this case, observe the contents of chapters "7.4 Cleaning surfaces" on page 46 and "7.2 Material compatibility" on page 45, as well as the chapter "2.2 Product-specific safety instructions" on page 32.
- ▶ If a local anaesthetic is being used, this should be thoroughly removed before the treatment.
- ▶ Always hold the handpiece firmly before switching it on or place it in the handpiece tray. If an unsecured handpiece is switched on, it may move uncontrollably through vibrations and prick or cause injury (see chapter "6.3 Switching the handpiece on and off" on page 43).

Injuries through contaminated needles or safety cartridges can result in the transmission of diseases (see chapter "5.7 Inserting or replacing the safety cartridge" on page 41):

- ▶ In the event of being injured by a contaminated needle, seek immediate medical attention from a physician!

In order to prevent the intensity of the treatment exceeding the extent desired:

- ▶ Prevent the risk of injury due to the needles protruding too deeply (see the "6.5 Setting the needle protrusion depth" on page 44 section).
- ▶ Avoid risk of injury due to the penetration frequency being too high (see chapter "6.4 Setting the penetration frequency" on page 44).

Please be sure to inform your patient to avoid contact with recently treated skin. Recently treated areas of skin should also be protected from

- Contamination and
- UV and solar radiation

. Additional irritants such as

- visits to swimming pools or saunas,
- abrasive or chemical peelings,
- hair removal at the treated areas or
- self-tanning products

should be avoided by the patient in the first two to three days after the treatment.

2.4 Contraindications and side effects

The following contraindications and side effects are the results of a thorough analysis of the professional clinical literature concerning microneedle treatment. If the person responsible for the treatment should have the smallest doubt that the safety of the patient cannot be guaranteed, for example due to secondary illnesses or conditions, the treatment should be withheld or stopped immediately.

2.4.1 Contraindications

With the following contraindications, **no** microneedling treatment should be performed:

- Haemophilia or other blood-clotting disorders
- Ongoing use of blood thinning medication (e.g. acetylsalicylic acid, warfarin, aspirin, heparin)
- Unregulated diabetes mellitus
- Any form of active acne in the area of the treatment
- Dermatitis (e.g. skin tumours, keloids or extreme susceptibility to keloid formation, solar keratosis, warts and/or moles) in the area of the treatment
- Open wounds and/or eczema and/or rashes in the area of the treatment
- Scars younger than 6 months
- Systemic infections and infectious diseases (e.g. hepatitis of types A, B, C, D, E or F; HIV infection) or acute local skin infections (e.g. herpes, rosacea)
- During chemotherapy, radiotherapy or high-dosage corticosteroid therapy (recommendation: from 4 weeks before the start until 4 weeks after the end of the therapy)
- Up to 12 months after a plastic surgery operation in the area to be treated
- Up to 6 months after filler injections in the area to be treated
- Allergy to topical anaesthetics (local anaesthetic)

- Contains nickel, not suitable for patients with known nickel allergies.
- May cause an allergic skin reaction.
- Under the influence of alcohol and/or drugs
- Pregnancy and lactation

The treatment of mucous membrane and eyeballs is strictly forbidden.

The treatment must be immediately aborted in the event of:

- Excessive perception of pain
- Fainting/dizziness

2.4.2 Side effects

In order to avoid unnecessary side effects:

- Find out the patient's medical history and possible previous surgical procedures before the treatment.
- Ensure that patients with a herpes simplex virus infection have been treated beforehand with a prophylactic antiviral therapy.

Complications during and after a microneedling treatment seldom occur in comparison to alternative methods of scar treatment and are generally minor. In most cases, this is localised bleeding in the area of the treated skin. In some cases, minor side effects may also arise. The following side effects are reported in the professional clinical literature.

Frequent:

- Pain and discomfort on the first day after the treatment
- Short-term inflammatory reactions, erythema and/or oedema up to 6 days after the treatment
- Skin irritations (e.g. itching or increased temperature), which normally subside during the first 12 to 72 hours after the treatment
- Formation of scabs, which normally recede in the first 5 days
- Temporary peeling of the skin which normally subsides within 8 days

Seldom:

- Formation of herpes simplex virus type I (HSV-I) blisters
- Formation of small pustules or milia as a result of inadequate cleaning of the skin prior to treatment
- Hyperpigmentation with the body's own pigments, in particular with darker skin types, but completely cleared up within a few weeks
- Inflammatory reactions, haematoma, erythema and oedema
- Retinoid reaction (from slight reddening to peeling of the skin)

As a matter of principle, recently treated skin areas should be protected from UV and solar radiation.

2.5 User requirements (qualifications)

The device and its accessories may only be used by persons who have acquired the following qualifications:

- Medical training
- Basic knowledge of microneedling treatment (see chapter "6.4 Setting the penetration frequency" on page 44 and "6.3 Switching the handpiece on and off" on page 43)
- Knowledge of health and safety provisions (see chapter "2.3 Important health and safety precautions" on page 33)
- Knowledge of risks and side effects (see chapter "2.4.1 Contraindications" on page 34 and "2.4.2 Side effects" on page 34)

2.6 Intended purpose, indications and proper use

2.6.1 Intended purpose

The intended purpose of the device is the percutaneous collagen induction (PCI) through minimally invasive injury to the epidermis and dermis by means of fine needle punctures (microneedling).

2.6.2 Indications

Treatment of facial acne scars for clinical improvement of the scarring, in particular through a reduction in the scar depth and the normalisation of the skin structure.

After a total of four treatments at intervals of 30 days respectively, it was possible to show a significant improvement in the acne scars on the basis of the Goodman and Baron Scar Scale within the scope of a clinical study performed by MT.DERM.

2.6.3 Proper intended use

The treatment must be carried out in a dry, clean and smoke-free environment as well as under hygienic conditions. The device must be prepared, used and looked after as described in this operating manual. In particular, the applicable requirements concerning workplace design must be observed and the materials to be used must be sterile.

Intended use also includes the assumption that this operating manual and in particular chapter "2 Important safety instructions" on page 32, have been read in full and understood.

Intended use also assumes the following consumer groups exclusively: Adult men and women who are in good health.

Improper use is when the device or its accessories are used in a manner other than as described in this operating manual or if the operating conditions have not been complied with. The treatment of patients with contraindications in particular is prohibited (see chapter "2.4.1 Contraindications" on page 34).

The device is designed to be placed on a level surface. During use, the handpiece is to be held in the user's hand in order to treat a sedentary or recumbent patient. In doing so, the user shall operate the device at a maximum distance of 70 cm.

2.7 Symbols on the product

The symbols described below can be found in this operating manual, on the device and on its accessories or the packaging:

SYM-BOL	MEANING
	Device complies with the requirements of the Medical Device Directive 93/42/EEG.
	Attention!
IP21	Housing offers protection from coarse dust and water drops
IPX6	Housing offers protection from jets of water
IP42	Housing offers protection from foreign objects thicker than 1 mm and against dripping water falling at an angle
	Applied part, type B: Applied part provides safety from electrical shock and leakage currents
	Handpiece
	Foot switch

SYMBOL	MEANING
	DC connection, inner pin positive
	DC
	AC
	Observe operating manual
	Manufacturer
	Date of manufacture
	Catalogue number, order number
	Serial number
	Lot code
	Sterilised with ethylene oxide
	Can be used until
	Temperature limit
	Humidity limit
	Air pressure limit
	Protect from moisture
	Only use indoors
	Fragile
	Do not use if packaging damaged

SYMBOL	MEANING
	Not for re-use
	Dispose of in compliance with waste electronic appliances guidelines
	Do not re-sterilise
	Warning of cutting or penetrating injuries
	Test plate: Inspection required in accordance with the manufacturer's instructions (safety check) every 24 months, see chapter "7.1 Routine maintenance" on page 45
	JAPAN TUV R-PSE
	Energy efficiency level IV
	Device complies with Japanese VCCI standards
	CHINA SJ/T 11364-2014
	Double insulation / device protection class II
RoHS 2	Device complies with the requirements of RoHS 2
	Ukraine UKRSeprO
	ETL mark – Intertek Testing Certification
	Device corresponds to the requirements of the S-Mark certificate EN 60950-1:2006 + A11:2009 + A1:2010 + A12:2011 + A2:2013

SYM-BOL	MEANING
	Device corresponds to the requirements of the Australia and New Zealand Regulatory Compliance
	Conformity of the power supply with FCC rule part 15, class B have been demonstrated with a standard output load.

SYM-BOL	MEANING
	Conformity of the power supply with Industry Canada, class B have been demonstrated with a standard output load.
	See operating manual for further information

3 SCOPE OF DELIVERY

- 1 Controller
- 1 Handpiece
- 1 Handpiece cable
- 1 Power supply
- 1 each of country-specific adapters (EU, UK, AU)
- 1 Foot switch
- 1 Handpiece tray

- 1 Operating manual
- 1 Declaration of Conformity
- 1 Inspection report for safety inspection ("STK report")

The safety cartridges and handpiece covers required for operation (E-0610) are not included in the delivery. Suitable accessory parts for the handpiece can be obtained via authorised retailers, see "4.3 Accessories and spare parts" on page 37.

4 PRODUCT INFORMATION

4.1 Technical data

Device type	AD61MP
Nominal voltage	15 V - DC
Power consumption	max. 27 VA
Power supply model	Model number: GTM96180-1817.9-2.9 Part number: WR9QG1200CSP-CR6B2958 (GlobTek)
Protection class	2
Penetration frequency	100-150 Hz (± 10%)
Drive	Precision motor - DC
Operating mode	Continuous operation
Dimensions (W x H x D)	196 mm x 52 mm x 175 mm
Handpiece weight	approx. 80 g
Total weight	approx. 1100 g
Applied part	Handpiece, type B
Needle protrusion depth	0 mm to 1.9 mm (± 0.15 mm)

4.2 Operating conditions

Ambient temperature	+15°C to +25°C
Relative humidity	30% to 75%
Air pressure	700 hPa to 1070 hPa

4.3 Accessories and spare parts

The following accessories and spare parts can be purchased from authorised retailers.

ACCESSORIES	ARTICLE DESCRIPTION	SCOPE OF DELIVERY	PIECE/PACKAGED UNIT
Safety cartridges	EMNPAED06MP (6-MNP 1.5-Exceed)	-	8
Handpiece cover for handpiece	E-0610	-	16

SPARE PART	ARTICLE DESCRIPTION	SCOPE OF DELIVERY	PIECE/ PACKAGED UNIT
Controller	AD61MP	1	1
Handpiece	CME21MP	1	1
Connecting handpiece cable	5E-G770	1	1
Foot switch	E1038MP	1	1
Handpiece tray	GAD60	1	1
Power supply	E116503	1	1

SPARE PART	ARTICLE DESCRIPTION	SCOPE OF DELIVERY	PIECE/ PACKAGED UNIT
Country-specific adapter	EU: E-1154-W	1	1
	UK: E-1155-W	1	1
	AU: E-1157-W	1	1
Operating manual	7EAD61MP	1	1
Declaration of Conformity	7EEGKAD61MP	1	1
Test Report for Inspection	–	1	1

5 GETTING STARTED

CAUTION

Cable tripping hazard

An incorrectly laid cable can constitute a tripping hazard and cause injury.

- ▶ Ensure all cables are positioned so that no one can trip over them or pull on them unintentionally.

CAUTION

Danger of injury and of device malfunctions

Unsuitable accessories and spare parts can impair the function and safety of the device. The device can be damaged or can fail or malfunction endangering persons.

- ▶ Use only the safety cartridges, accessories and spare parts listed in chapter "4.3 Accessories and spare parts" on page 37.

CAUTION

Risk of short circuit and risk of light electric shocks

There is a risk of electrical shocks and damage to electronics if there is visible damage to cables or cable connections.

- ▶ Check the device and the cables by visually checking for damage such as a defective cable connection for example.
- ▶ Never kink the product's cable.

ATTENTION

Damage from condensation

If the device is exposed to high temperature variations e.g. during transport, condensate may accumulate inside and damage the electronics.

- ▶ Ensure that the device has reached ambient temperature before using it. If the device has been exposed to high temperature variations, wait for at least 3 hours per 10 °C of temperature difference before using it.
- ▶ Use the device only with an ambient temperature of +15°C to +25°C.

5.1 Positioning the device

⚠ CAUTION

Reduced functionality due to electromagnetic interference

Portable and mobile HF communications devices, such as mobile phones or WLAN routers, can influence the functionality of the device through the transmission of electromagnetic radiation. Safe operation of the device can no longer be assured.

- ▶ Avoid sources of interference by ensuring that portable or mobile radio devices are not operated in the vicinity of the device.
- ▶ Inform your patients about the risk of sources of interference.

ATTENTION

Damage to the device through inadequate stability

If the device is not positioned safely, safe operation is not guaranteed. The device could fall and become damaged during operation.

- ▶ Place the device on a clean, firm and level surface.
 - ▶ Ensure that the buttons and the display of the device remain freely accessible during use.
 - ▶ Never place the device on or underneath other devices.
-
- ▶ Position the device such that the recommended safety distance from possible sources of interference is maintained. You can find information on the recommended safety distances in the chapter Information on electromagnetic compatibility.
 - ▶ Position the device so that the power supply is accessible at all times.

5.2 Connecting the power supply

The device may only be operated with the power supply cited in the operating manual of the device. The mains voltage must match with the device voltage cited on the label of the power supply.

- ▶ If there is no suitable country-specific adapter included in the scope of delivery, contact an authorised dealer (see also chapter "10 Questions and problems" on page 47).

To remove the existing country-specific adapter, if necessary:

- ▶ Pull back and hold the catch for the adapter and at the same time lift the country-specific adapter out of the recess in the power supply.


Fig. 1: Remove country-specific adapter

To insert the desired country-specific adapter:

- ▶ Place the country-specific adapter in the recess of the power supply as shown below.
- ▶ Press the adapter into the plug-in power unit until the catch can be heard to engage.


Fig. 2: Inserting country-specific adapter

In order to establish the power supply:

- ▶ Insert the DC plug of the power supply into the socket for the power supply (⊕ ⊖) on the rear of the controller.


Fig. 3: Socket for the power supply

- ▶ Insert the power supply into a mains outlet. The standby LED illuminates blue.
- ▶ Make sure that the power supply is accessible at all times.

5.3 Connecting the foot switch (optional)

Only the foot switch cited in chapter "4.3 Accessories and spare parts" on page 37 will ensure safe operation.

At first delivery, the "CONTINUAL" footswitch mode is selected (see chapter "6.2.2 Selecting the foot switch mode" on page 43).

- ▶ Insert the jack plug for the foot switch into the socket for the foot switch (↘) on the back of the control unit.


Fig. 4: Socket for the foot switch

5.4 Connecting the handpiece

If no handpiece is connected to the controller, a message appears on the display:

Fig. 5: Handpiece detached

- ▶ Connect the handpiece cable jack socket supplied in the scope of delivery to the connection on the handpiece of the device.
- ▶ Insert the FireWire plug of the handpiece cable fully into the output socket of the controller.


Fig. 6: Output socket

The handpiece is detected by the controller and the pre-set puncturing frequency appears on the display.

If the handpiece is not detected:

- ▶ Refer to chapter "10 Questions and problems" on page 47.

5.5 Cleaning and disinfecting the equipment

ATTENTION

Damage due to fluids

If cleaning products or disinfectant gets inside the handpiece, the controller or the power supply, this may cause a short-circuit. In addition, this may also corrode the electrical and mechanical components.

- ▶ Never immerse the handpiece, the controller or the power supply in cleaning product or disinfectant.
- ▶ Never clean the handpiece, the controller or the power supply in autoclaves or in an ultrasonic cleaning bath.
- ▶ Consider material compatibilities when selecting cleaning materials and disinfectants (see chapter "7.2 Material compatibility" on page 45).

ATTENTION

Device damage caused by unapproved cleaning materials or disinfectants

Cleaning materials or disinfectants that are not compatible with the materials of the device can cause damage to its surface.

- ▶ Use only cleaning materials and disinfectants that are authorised for use by the legislation in your country.
- ▶ Consider material compatibilities when selecting cleaning materials and disinfectants (see chapter "7.2 Material compatibility" on page 45).
- ▶ Disconnect the controller from the power supply.
- ▶ Wipe the controller, the handpiece cable, the handpiece and the handpiece tray with a soft cloth moistened with cleaning product.
- ▶ Wipe the handpiece, the controller, the handpiece cable and the handpiece tray with a soft cloth moistened with disinfectant.

5.6 Pulling the handpiece cover on

We recommend using the handpiece covers provided as accessories. Alternatively, you can use equivalent handpiece covers with a diameter of 20 to 27 mm. This represents a sleeve width of approx. 30 to 45 mm.

In order to protect the handpiece and the handpiece cable from contamination with bodily fluids or accompanying preparations contaminated with bodily fluids during the treatment:

- ▶ Unplug the handpiece cable from the handpiece.
- ▶ Slide the handpiece cover fully over the handpiece cable.
- ▶ Connect the handpiece cable, with the handpiece cover fitted on it, to the handpiece connection once more.
- ▶ Draw the end of the handpiece cover fully over the handpiece.
- ▶ Stick the end of the handpiece cover firmly to the moving part of the handpiece.


Fig. 7: Handpiece cover in place

- ▶ Dispose of the used hygiene sleeves in containers for proper disposal of contaminated parts, after the conclusion of each treatment.

5.7 Inserting or replacing the safety cartridge

WARNING

Danger of injury due to needles

Injuries through contaminated needles can result in the transmission of diseases.

- ▶ Only use original amiea med safety cartridges with the device that match the supplied handpiece. The safety membrane integrated in the safety cartridges prevents the return flow of (body) fluids as well as the ingress of pathogens into the handpiece.
- ▶ Never touch the tip of the safety cartridge.
- ▶ Switch the handpiece off before replacing the safety cartridge.
- ▶ Never press the needles out of the safety cartridge.
- ▶ When removing a safety cartridge, ensure that any fluids such as ink or accompanying preparations that may be being used do not run from the safety cartridge into the handpiece.
- ▶ After removing the used safety cartridge, check that the needles have been completely withdrawn into the cartridge.
- ▶ Dispose of used safety cartridges that are no longer required for the current treatment properly (see chapter "9 Disposal" on page 46).
- ▶ The safety cartridges are shipped in sterile packaging. Only remove them from the sterile packaging immediately before use. To do so, open the sterile packaging of the safety cartridge at the PEEL marking ().
- ▶ Always place the handpiece back in the handpiece tray so that the safety cartridge does not rest on the surface.
- ▶ In the event of being injured by a contaminated needle, seek medical attention from a physician.

The needles of the amiea med safety cartridges are produced from stainless steel (1.4310, X10CrNi18-8, AISI 301). The tips of the amiea med safety cartridges are produced from PC (polycarbonate).

- ▶ Switch the handpiece off.
- ▶ Hold the safety cartridge by the base body and when inserting it in the handpiece (never hold it by the tip!).

- ▶ Insert the safety cartridge into the threaded connection of the handpiece. You can insert the safety cartridge with all possible orientations.


Fig. 8: Inserting the safety cartridge

- ▶ Turn the safety cartridge clockwise until you can feel resistance. The safety cartridge is now firmly seated in the handpiece.


Fig. 9: Tightening the safety cartridge

- ▶ Check that the safety cartridge is securely attached by pulling gently on it.
- ▶ To remove the safety cartridge, turn it anticlockwise and remove it from the handpiece.

5.8 Checking equipment

- ▶ Carry out a visual check of the equipment:
 - Are there any signs of external damage (e.g. kinks in the handpiece cable, detached cap in the area of the cable connection)?
 - Are any cables disconnected?
 - Are the safety cartridge and needles aligned correctly? Are all needles fully withdrawn into the safety cartridge ?
 - Has the date for the next inspection (biennial service) on the test plate of the controller passed?
 - Does the controller display an error message?
- ▶ Switch the handpiece on and carry out an acoustic check: Are operating noises or sounds prominent?
- ▶ In the event of anomalies, refer to chapters "10 Questions and problems" on page 47 and "11.1 Guarantee" on page 48.
- ▶ If the device can no longer be operated harmlessly, take it out of service and contact an authorised dealer.

6 USE

6.1 Switching the control unit to operating mode or standby mode

To switch the controller to operating mode:

- ▶ Ensure that the power supply is established (see chapter "5.2 Connecting the power supply" on page 39).
- ▶ Press the **Standby/On** button (⏻). The standby button illuminates.

The start screen appears on the display.

If a handpiece is connected, the default puncture frequency (120) appears.

If an error message appears:

- ▶ Refer to the chapter "10 Questions and problems" on page 47.


To switch the controller to standby mode:

- ▶ Press the **Standby/On** button (⏻) again. The display goes out and the standby LED illuminates blue.

6.2 Configuring the controller (menu)

To enter the menu:

- ▶ Simultaneously press both the **left** (◀) and **right** (▶) buttons on the controller. The "KEY TONES" menu item appears on the display.

To toggle between the three menu items:

- ▶ Press the **right** button (➤) to navigate to the next menu item.
- ▶ Press the **left** button (➤) to navigate to the previous menu item.

To exit the menu:

- ▶ Switch to the "EXIT" menu item.

- ▶ Press the **SET** button (⌚).

Alternative:

- ▶ Press the **Standby/On** button (⌚).

Alternative:

- ▶ Simultaneously press both the **left** (➤) and **right** (➤) buttons.


Alternative:

- ▶ Wait approx. one minute until the menu automatically closes.

6.2.1 Switching the key tones on and off

You can switch the key tones on and off.

- ▶ Switch to the "KEY TONES" menu item.


- ▶ Press the **SET** button (⌚) to change the setting for the key tones ("ON" or "OFF").

6.2.2 Selecting the foot switch mode

You can select between two modes for operating the foot switch:

- "ON/OFF" (switch mode): The handpiece is alternately switched on and off again when pressing the foot switch.
- "CONTINUAL" (jog mode): The handpiece remains operating until the foot switch is pressed.

- ▶ Switch to the "FOOT SWITCH MODE" menu item.

- ▶ Press the **SET** button (⌚) to change the setting for the foot switch mode.

6.3 Switching the handpiece on and off

CAUTION

Danger of injury when using the handpiece

If an unsecured handpiece is switched on, it may move uncontrollably through vibrations and cause injury.

- ▶ Hold the handpiece firmly before switching it on or place it in the handpiece tray.

To switch on the handpiece:

- ▶ Ensure that the power supply is established (see chapter "5.2 Connecting the power supply" on page 39).
- ▶ Ensure that the controller is switched to operating mode (see chapter "6.1 Switching the control unit to operating mode or standby mode" on page 42)
- ▶ Ensure that the handpiece is correctly connected (see chapter "5.4 Connecting the handpiece" on page 40).
- ▶ Press the **SET** button (⌚) on the controller. Ongoing handpiece operation is indicated on the controller display by a rotating symbol (⌚).

If the connected handpiece does not start:

- ▶ Observe the chapter "10 Questions and problems" on page 47.

To switch off the handpiece:

- ▶ Press the **SET** button (⌚) again.

Alternatively, you can switch the handpiece on and off with the foot switch. The foot switch function depends on the foot switch mode setting in the device menu (see chapter "6.2.2 Selecting the foot switch mode" on page 43).

When using the switch mode ("ON/OFF"):

- ▶ Briefly press the foot switch to switch the handpiece on or off.

When using the jog mode ("CONTINUAL"):

- ▶ Only step on the foot switch as long as you want to operate the handpiece.

6.4 Setting the penetration frequency

⚠ CAUTION

Danger of injury due to the penetration frequency being too high

If the penetration frequency is too high the skin can be cut

- ▶ Always start with a low puncture frequency, eg. 100 punctures per second.
- ▶ Select the penetration frequency on the basis of the patient's skin characteristics and suitable for the working speed.

The preset penetration frequency is 120 penetrations per second. The controller automatically resets to this preset penetration frequency every time it is switched on. The penetration frequency can be 11 stages of 100-150 penetrations per second.

The number of penetrations per second shown on the display may deviate by $\pm 10\%$ in stationary operation from the value that actually occurs at the handpiece.

To set a suitable penetration frequency:

- ▶ Ensure that the controller is switched to operating mode (see chapter "6.1 Switching the control unit to operating mode or standby mode" on page 42)
- ▶ Press the **right** button (➤) to increase the penetration frequency.
- ▶ Press the **left** button (◀) to decrease the penetration frequency.

6.5 Setting the needle protrusion depth

⚠ CAUTION

Danger of injury due to excessive needle protrusion depth

Greater needle protrusion depth allows a greater penetration depth. If the needle protrusion depth and penetration depth are too great, the risk of side effects, such as haematoma, erythema or oedema, increases (see chapter "2.4.2 Side effects" on page 34).

- ▶ Choose the needle protrusion depth depending on the skin condition of the patient, the area of the skin to be treated and the pressure applied.
- ▶ The needle protrusion depth shown on the scale may deviate by ± 0.15 mm from the actual value.

The moving part of the handpiece is labelled with a scale, which indicates the needle protrusion depth setting. The needle protrusion depth can be adjusted between 0 mm and 1.9 mm.

- ▶ Switch the controller to operating mode.
- ▶ Switch the handpiece on.
- ▶ Set the needle protrusion depth by turning the moving part of the handpiece (1). Use the scale as a guide and turn it until the marker (2) is at the desired position. Turning it clockwise (when looking at the needles) increases the needle protrusion depth.

Clinical studies by MT.DERM have shown that the decisive factor in effective acne scar treatment is the induction of so-called pinpoint bleeding. Treatment should be started with the smallest possible needle protrusion depth with the depth being successively increased until pinpoint bleeding occurs. As soon as point bleeding occurs in the treated area, treatment should be switched to the next area.


Fig. 10: Setting the needle protrusion depth

7 CLEANING AND MAINTENANCE

CAUTION

Risk of short circuit and risk of light electric shocks

There is a risk of minor electrical shocks and damage to electronics on the power unit if device components are dismantled and cleaned whilst still electrically connected.

- ▶ Disconnect the controller from the power supply by removing the mains adapter from the socket before any cleaning or maintenance tasks.

ATTENTION

Damage due to fluids

If cleaning products or disinfectant gets inside the handpiece, the controller or the power supply, this may cause a short-circuit. In addition, this may also corrode the electrical and mechanical components.

- ▶ Never immerse the handpiece, the controller or the power supply in cleaning product or disinfectant.
- ▶ Never clean the handpiece, the controller or the power supply in autoclaves or in an ultrasonic cleaning bath.
- ▶ Consider material compatibilities when selecting cleaning materials and disinfectants (see chapter "7.2 Material compatibility" on page 45).

7.1 Routine maintenance

As part of the routine maintenance of the device, an inspection (safety inspection) is required every 24 months. Inspections may only be carried out by authorised dealers and are to be performed according to manufacturer instructions. All necessary inspection steps for the inspection of the device are stated in the inspection report from the scope of delivery.

Return the device to a specialist retailer for inspection in accordance with the manufacturer's instructions (safety check), at regular 24 month intervals.

- ▶ Before use, always ensure that the date for the next inspection on the test plate of the control unit has not passed.
- ▶ Do not use the device if the date on the test plate has passed and commission a new inspection if necessary.

7.2 Material compatibility

- ▶ To clean the device, use a mild soap detergent authorised in your country, or 50% diluted 1-propanol solution.
- ▶ Use the following disinfectants to disinfect the device:

MANUFACTURER	PRODUCT	WORKING-IN-TIME
Antiseptica	Kombi Liquid	1 to 5 min
Bode Chemie	Bacillol	30 s to 1 min
Ecolab	Incidin Foam	1 to 2 min
Schülke & Mayr	Mikrozid Liquid	1 to 2 min

Based on the material compatibilities, you can assess which cleaning products and disinfectants are suitable.

The device surfaces are **resistant** to:

- Mild acids (e.g. boric acid $\leq 10\%$, acetic acid $\leq 10\%$, citric acid $\leq 10\%$)
- Aliphatic hydrocarbons (e.g. pentane, hexane)
- Ethanol
- Most inorganic salts and their aqueous solutions (e.g. sodium chloride, calcium chloride, magnesium sulphate)

The device surfaces are **not resistant** to:

- Strong acids (e.g. hydrochloric acid $\geq 20\%$, sulphuric acid $\geq 50\%$, nitric acid $\geq 15\%$)
- Oxidising acids (e.g. peracetic acids)
- Lyes (e.g. caustic soda, ammonia and all compounds with $\text{pH} > 7$)
- Aromatic/halogenated hydrocarbons (e.g. phenol, chloroform)
- Acetone and petrol

7.3 Disinfecting surfaces

Before and after each use:

- ▶ Remove the safety cartridge from the handpiece and pull the handpiece cable out of the unit.
- ▶ Disinfect the handpiece cable, the handpiece and the handpiece tray as described in chapter "5.5 Cleaning and disinfecting the equipment" on page 40.

7.4 Cleaning surfaces

External soiling:

- ▶ Remove the safety cartridge from the handpiece and pull the handpiece cable out of the unit.
- ▶ Wipe the device and its accessories with a soft cloth moistened in the cleaning product or disinfectant.

8 TRANSPORT AND STORAGE CONDITIONS

8.1 Controller, handpiece with handpiece cable, power supply and foot switch

ATTENTION

Damage due to dropping

If the controller or the handpiece is dropped it can be damaged.

- ▶ Always set the controller down on the rubber feet on a flat and stable surface so that it cannot fall.
- ▶ Place the handpiece back in the handpiece tray after each use so that it is safe.
- ▶ If the controller or the handpiece is dropped, carry out a visual check of the components.
- ▶ Hand the controller or the handpiece in to an authorised dealer if the components have visible damage or if they no longer function normally.

- ▶ Transport the device and its accessories only in the original packaging.
- ▶ Store the device and its accessories only in the following conditions:

Ambient temperature	-10°C to +40°C
Relative humidity	10% to 85%
Air pressure	540 hPa to 1070 hPa

8.2 Safety cartridge

- ▶ Always transport and store the safety cartridges in the following conditions:

Ambient temperature	+15°C to +25°C
Relative humidity	30% to 65%

9 DISPOSAL

The product service life foreseen by the manufacturer encompasses 5 years for all parts included in the delivery insofar as the operating, storage and transport conditions described in this operating manual as well as all maintenance tasks are complied with.

If the product and its accessories continue to be used after this product service life has expired, there is a risk of increased wear or the failure of individual functions. Fault-free operation will no longer be assured.

After the product service life has expired:

- ▶  Dispose of the product showing the symbol opposite, in compliance with applicable regulations on the disposal of electronic waste. Make enquiries with your specialist retailer or authorities about the applicable regulations.
- ▶ Place a container at your workplace for the proper disposal of contaminated parts.
- ▶ Dispose of used or faulty safety cartridges and safety cartridges, whose packaging is damaged in a non-penetrable sharps container in accordance with the regulations applicable in your country.

10 QUESTIONS AND PROBLEMS

The following functional disruptions or fault messages could occur:

FAULT/MESSAGE	FAULT	CAUSE	RECTIFICATION
The controller will not switch to standby mode.	The power supply is disconnected.	The power supply is not correctly connected.	<ul style="list-style-type: none"> ▶ Check the connection of the power supply to the controller (see chapter "5.2 Connecting the power supply" on page 39). ▶ Check the power supply of the work station.
The controller switches off automatically during operation.	The power supply is disconnected or there is a fault in the controller.	The power supply is not correctly connected.	<ul style="list-style-type: none"> ▶ Check the connection of the power supply to the controller (see chapter "5.2 Connecting the power supply" on page 39). ▶ Check the power supply of the work station. <p>If the power supply is connected:</p> <ul style="list-style-type: none"> ▶ Switch the controller to standby mode and wait for 30 seconds (see chapter "6.1 Switching the control unit to operating mode or standby mode" on page 42).

FAULT/MESSAGE	FAULT	CAUSE	RECTIFICATION
"Handpiece or device overloaded!"	There is an error in the handpiece.	The handpiece cable is not connected correctly.	<ul style="list-style-type: none"> ▶ Check the connections between the handpiece and controller. ▶ Press any button on the controller to return to the default display. ▶ Check that the safety cartridge is securely attached to the handpiece. ▶ Press any button on the controller to return to the default display.
"Handpiece detached!"	The handpiece is not detected.	The handpiece cable is not connected correctly	<ul style="list-style-type: none"> ▶ Carry out all of the steps in the chapter "5.4 Connecting the handpiece" on page 40.
"System Error!"	There is an error in the handpiece or controller.		<ul style="list-style-type: none"> ▶ Switch the controller to standby mode and wait for 30 seconds (see chapter "6.1 Switching the control unit to operating mode

FAULT/MESSAGE	FAULT	CAUSE	RECTIFICATION
			<p>or standby mode" on page 42).</p> <ul style="list-style-type: none"> ▶ Switch the controller back to operating mode. ▶ If it is not possible to establish standby mode, disconnect the power supply

FAULT/MESSAGE	FAULT	CAUSE	RECTIFICATION
			<p>from the controller for 30 seconds.</p> <p>In order to rectify further functional disruptions:</p> <ul style="list-style-type: none"> ▶ Firstly, disconnect all components, such as the handpiece and the power supply, from the controller. ▶ Check all connections and reconnect the components. ▶ Check the function of the device again. ▶ If malfunctions persist or if you have any questions or complaints, please contact your specialist retailer. <p>For our current offers, a selection of safety cartridges and accessories, visit www.amieamed.com.</p> <p>Visit www.amieamed.com/contact-us/distributors to find your local dealer and related contact information.</p>

11 MANUFACTURER DECLARATIONS

11.1 Guarantee

This device is a high-quality medical device.

The product is covered by the statutory 2-year warranty against malfunctions due to material failure or production defects.

We guarantee that the safety cartridges will remain sterile for 5 years in sealed and undamaged packaging, if the transport and storage conditions are complied with. In the event of complaints relating to safety cartridges, please provide us with the lot number printed on the label.

The warranty does not apply to the following damage:

- Damage and consequential damage caused by improper use or failure to comply with the operating manual.
- Damage caused by the ingress of fluids or dirt into the inside of the handpiece, the controller or the power supply.

11.2 Electromagnetic compatibility

WARNING

The device has been tested for electromagnetic compatibility (EMC) and is suitable for use in the area of domestic healthcare (CISPR 11 class B). Use in motor vehicles and aeroplanes is not permitted.

Caution when using this device next to or in connection with other devices. This could lead to defective operation.

If such use or installation is necessary, this device and the other devices must be checked to ensure they are functioning normally.

When using accessory parts, power units and cables that are not included in the "4.3 Accessories and spare parts" on page 37 of the device manufacturer, there is a risk of increased electromagnetic emissions or reduced electromagnetic immunity of this device and therefore defective operation.

Portable HF communication devices (such as telephones, tablets, etc.) should not be used at distances of less than 30 cm from the device and its electrical accessories. Otherwise this could impair/adversely affect the capacity of the device.

With regards to guaranteeing its basic safety and significant performance characteristics, the device has been tested for its electromagnetic immunity and its electromagnetic radiation. These remain within the limit values specified in the following tables, with consideration to the stipulated minimum distance from neighbouring devices. No maintenance measures are required to guarantee basic safety and the significant performance characteristics.

11.2.1 Manufacturer declarations on electromagnetic emissions

The device is intended for operation in the electromagnetic environment specified below. The user of the device should ensure that it is operated in this type of environment.

EMISSIONS MEASUREMENTS	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDELINES
HF emissions per CISPR 11	Group 1	The device uses HF energy exclusively for its internal function. Its HF emissions are therefore very low and are not likely to cause any interference in nearby electronic equipment.
HF emissions per CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to a public power supply network that supplies buildings also used for domestic purposes.
Harmonics per IEC 61000-3-2	not applicable	
Voltage fluctuations/flicker per IEC 61000-3-3	not applicable	

11.2.2 Manufacturer declarations on electromagnetic interference immunity

The device is intended for operation in the electromagnetic environment specified below. The user of the device should ensure that it is operated in this type of environment.

IMMUNITY TESTS	IEC 60601-TEST LEVEL	COMPLIANCE LEVEL	ELECTRO-MAGNETIC ENVIRONMENT - GUIDELINES
Discharge of static electricity (ESD) per IEC 61000-4-2	±8 kV contact discharge ± 2, 4, 8, 15 kV air discharge	in compliance	Floor should be covered with ceramic tiles. If the floor is covered with synthetic material, the relative air humidity must be at least 30%.
Rapid transient electrical dis-	±2 kV for mains cables	in compliance	The quality of the supply voltage should

IMMUNITY TESTS	IEC 60601-TEST LEVEL	COMPLIANCE LEVEL	ELECTRO-MAGNETIC ENVIRONMENT - GUIDELINES
turbances/bursts per IEC 61000-4-5	± 1 kV for incoming and outgoing lines		reflect a typical business or hospital environment.
Impulse voltage (surges) per IEC 61000-4-5	± 1 kV differential mode voltage ±2 kV common-mode voltage	in compliance	The quality of the supply voltage should reflect a typical business or hospital environment.
Voltage dips, brief interruptions and fluctuations in the supply voltage per IEC 61000-4-11	0% UT; 0.5 periods at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 periods and 70% UT 25/30 periods. Single phase: at 0°	in compliance	The quality of the supply voltage should reflect a typical business or hospital environment. If the user of the device requires continued operation also with interruptions to the power supply, it is recommended that the device be supplied by an uninterruptible power supply or a battery.
Voltage interruption	0% UT, 250/300 periods *UT is the mains alternating voltage before the application of the test level.	in compliance	
Magnetic field with the supply frequency (50/60 Hz) per IEC 61000-4-8	30 A/m at 50 Hz and at 60 Hz	in compliance	Magnetic fields with the mains frequency should reflect the typical values, as found in a business or hospital environment.

IMMUNITY TESTS	IEC 60601-TEST LEVEL	COMPLIANCE LEVEL	ELECTRO-MAGNETIC ENVIRONMENT - GUIDELINES
----------------	----------------------	------------------	---

Immunity in compliance with conducted interference. IEC 61000-4-6

3 V_{eff}
150 kHz to 80 MHz
6 V_{eff} within the ISM band and the amateur frequency bands between 0.15 MHz and 80 MHz

Portable and mobile radio devices are not used at a shorter distance from the device including the cables, than the recommended protective distance, which is calculated in accordance with the equation suitable for the transmission frequency. With P as the rated power of the transmitter in Watt (W) in accordance with the specifications of the transmitter manufacturer and d as the recommended protective distance in metres (m). With 80 MHz and 800 MHz the higher value applies. These guidelines may not apply in all situations. The propagation of electromagnetic waves is influenced by their absorption by and reflection

IMMUNITY TESTS	IEC 60601-TEST LEVEL	COMPLIANCE LEVEL	ELECTRO-MAGNETIC ENVIRONMENT - GUIDELINES
----------------	----------------------	------------------	---

Immunity against high frequency electromagnetic fields IEC 61000-4-3

10 V/m, 80MHz to 2.7 GHz

in compliance

off buildings, objects and people.

The field strength of static radio transmitters determined through local analysis* should be lower at all frequencies than the compliance level**. Interference is possible in the vicinity of devices that bear the following symbol.

*The field strength of static transmitters, such as base stations of mobile phones and land mobile radio services, amateur stations, AM and FM radio and television broadcasting stations cannot be theoretically precisely predetermined. In order to determine the electromagnetic environment with regards to static transmitters, a study of the location should be considered. If the determined field strength at the location of the device exceeds the compliance level specified above then the device must be observed in its normal operation, to verify the intended function. If unusual performance characteristics are observed, additional measures may be necessary, such as a change in the orientation or in the device location.

**Over the frequency range from 150 kHz to 80 MHz, the field strength is less than 3 V/m.

ISM band (Industrial, Scientific and Medical Band)

FREQUENCY [MHZ]	FREQUENCY BAND [MHZ]	MODULATION	IMMUNITY TESTING [V]
6.78	6.765 to 6.795	AM 80% 1 kHz	6
13.56	13.553 to 13.567	AM 80% 1 kHz	6
27.12	26.957 to 27.283	AM 80% 1 kHz	6
40.68	40.66 to 40.7	AM 80% 1 kHz	6

Amateur radio

FREQUENCY BAND [MHZ]	MODULATION	IMMUNITY TESTING (V)
1.8 to 2.0	AM 80% 1 kHz	6
3.5 to 4.0	AM 80% 1 kHz	6
5.3 to 5.4	AM 80% 1 kHz	6
7 to 7.3	AM 80% 1 kHz	6
10.1 to 10.15	AM 80% 1 kHz	6
14 to 14.2	AM 80% 1 kHz	6
18.07 to 18.17	AM 80% 1 kHz	6
21.0 to 21.4	AM 80% 1 kHz	6
24.89 to 24.99	AM 80% 1 kHz	6
28.0 to 29.7	AM 80% 1 kHz	6
50.0 to 54.0	AM 80% 1 kHz	6

Test requirements for immunity against wireless radio equipment

TEST FREQUENCY	FREQUENCY BAND [MHZ]*	RADIO SERVICE*	MODULATION**	MAX. POWER [W]	DIS-TANCE [M]	IM-MUNI-TY TESTING [V/M]
385	380 to 390	TETRA400	Pulse modulation** 18 Hz	1.8	0.3 (12 inch)	27
450	430 to 470	GMRS 460, FRS 460	FM*** ± 5 kHz deviation, 1 kHz Sinus	2	0.3 (12 inch)	28
710	704 to 780	LTE band 13, 17	Pulse modulation** 217 Hz	0.2	0.3 (12 inch)	9
810	800 to 930	GSM 800/900, TETRA A 800, iDEN820, CDMA 850, LTE band 5	Pulse modulation** 217 Hz	2	0.3 (12 inch)	28
1720	1700 to 1970	GSM 1800, CDMA 1900, GSM 1900, DECT	Pulse modulation** 217 Hz	2	0.3 (12 inch)	28

TEST FREQUENCY	FREQUENCY BAND [MHZ]*	RADIO SERVICE*	MODULATION**	MAX. POWER [W]	DISTANCE [M]	IMMUNITY TESTING [V/M]
----------------	-----------------------	----------------	--------------	----------------	--------------	------------------------

LTE
band
1, 3, 4,
25,
UMTS

2450	2400 to 2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE band 7	Pulse modulation** 217 Hz	2	0.3 (12 inch)	28
------	--------------	---	---------------------------	---	---------------	----

5240	5100 to 5800	WLAN 802.11 a/n	Pulse modulation** 217 Hz	0.2	0.3 (12 inch)	9
5500	5785					

* With some radio services, only the Uplink channel is considered.

** The carrier is modulated with a square wave signal with 50% duty cycle.

*** As an alternative to FM modulation, pulse modulation can be used of 50% at 18 Hz, because this would be the worst case although this does not constitute actual modulation.

The device must be set up and put into operation with consideration to the EMC information, this applies in particular to the distances from other devices emitting radio waves. Never use the device at a distance of less than 30 cm from another device emitting radio waves.

In the event of unlikely interference, which could be caused by devices emitting radio waves, observe the minimum distances and move the device away from any radio waves if necessary.

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MT.DERM GmbH | Blohmstraße 37-61 | 12307 Berlin / Germany